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This series of Lifumat furnaces has been developed to offer many advantages to Metallurgical Research and Quality Control environments within industry when compared with conventional techniques. These casting machines are very versatile and capable of a wide range of melting applications from platinum to scrap iron with excellent

reproducibility, minimal elemental losses, high sample throughput, ease of use and reliability. As many as 250 samples per day may be prepared with 24 hour operation. Their high melting capacity with low energy consumption, reproducible operation through microprocessor control, ease of maintenance and a limited space requirement are the most outstanding features of this series. Lifumat-Met-6.6, -7.7, Super and their vacuum versions are also suitable for production of calibration standards (cross sectional samples) and scrap (recycling) up to 1 kg.

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


Vakuumschleuderarm / Vacuum centrifugal arm 3.3 Vac

Remelting with a LINN unit by induction heating provides the metallurgical laboratory in research, production and recycling industries with a system for single and routine testing which offers many advantages compared with traditional procedures of sample preparation. The important advantages are **reproducibility**, **low element losses**, as well as **high throughput** which are confirmed by the installation of many units for all kinds of application from precious metals to scrap recycling.

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The Lifumat-1.2-Table Ox offers rapid production of spectrometer samples by induction heating. Ideal for sample preparation of oxydic materials for XRF/RFA, AAS and ICP analytical systems. This compact unit will satisfy the fusion needs of small laboratories. The 10 step selectable melting power is generated by a powerful solid state power generator. A crucible vibration facility may be operated as required. This small table top unit is simple to install and safe to operate with low power consumption.

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

The well established fusion system, the Lifumat-2.0 Ox is available with a 25 program controller / programmer. A tilting and turning device improves the homogeneity and can also be used for high melting materials. Fusion temperatures up to 1500 °C will be achieved. An adjustable cooling system ensures best cooling of the sample.

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

High-frequency induction heating

By this method of heating beads fused in Pt/Au or graphite crucibles with temperatures up to 1500 °C can be achieved within 3 min. The advantage of **induction fusion** is easy handling and **operation safety** at increased temperature with **minimum time**. The eddy currents guarantee that the melt has an optimal mix by the so-called bath movement. By that a constant and homogeneous quality of the melt is achieved. Accuracies of analysis and economy with optimal sample quantity confirm the high quality of this method with a minimum of elemental losses, which can be rectified. Graphite-, glass-carbon-, Zirconia- and Pt/Au-crucibles can be applied. Other fusion methods as i.e. heating by torch or within resistance-heated furnaces show many disadvantages: unequal temperature distribution, poor mixing, environmental problems, high operation costs and larger fusions.

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Semi-automatic high frequency fusion furnace for spectrometer sample preparation

This induction heated fusion system allows the laboratory to produce quickly and economically up to 35 reproducible fusions of oxidic materials in one step. The melting conditions can be selected individually for each sample. This provides a basis for a constant quality in the following evaluation of the fusion for XRF/RFA, AAS, ICP. The main requirements to fusion are: reproducibility, time saving, costs effective, environmental acceptability, no gas, flexibility, safe operation, high temperature, ergonomics. These are fulfilled to a high degree by this system. Special attention was given to a high sample throughput and melting power at low energy consumption, simple operation and Simatic-SPS-control and easy servicing by board - replacement, low space requirement. Also, the crucible materials are not limited to Pt/Au. Graphite, glass-carbon or met. zirconia can also be used. Typical applications are: cement industry, mines, foundries etc. The system can be equipped with a separate PC for complete control and logging of processed samples. The Windows 95-based control program visualizes the unit and entry of sample preparation parameters and sample identification is intuitive and easy. All processed samples are documented in a file with achieved temperature if a pyrometer is installed.

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Induction Fusion and Remelt Furnace

This new sample preparation unit enables laboratories to carry out both methods of sample preparation in one unit - **remelting of metals and induction fusion of oxidic samples for XRF/RFA, AAS, ICP, X-Emission**.

The Lifumat-C-2000-3.3 works with a powerful, safe to operate high frequency generator, coupled with a high frequency power switchover for both methods. Control is done by a microprocessor control with interface for a master computer.

The sample preparation unit can be switched over between manual and automatic operation. Measurement and temperature control by optical pyrometer.

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max. weight of casting and samples, sample sizes:

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